


## Cassette tray 126

X-ray cassette tray for different X-ray cassette sizes. Slides on the X-ray cassette rails under the table. Stainless steel.

- Size 395 x 595 mm
- Stainless steel
- Suitable with Promerix and Practico



### Technical Specifications

Weight	2.2 kg
Product Code	100000126

### Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)