

Multipurpose accessory cart 60220

Accessory cart for mobile storage and moving Merivaara operating table sections and accessories. The cart helps to keep operating room more organized and reduces set up times.

- Size 580 x 810 mm, height 1770 mm
- 4 castors with lockable brakes
- Holds up to total of 12 sections on both sides of the cart
- Includes 6 accessory rails which enable variable accessory storing
- 2 baskets 590 x 300 x 150 mm can be placed in 6 different locations
- Long sections, such as lightweight transfer board and carbon fibre backrests, can be placed in vertical position in two locations
- Narrow head rests can be stored using 100 060 146 adaptor is not included
- Accessory rails available: EUR, US, UK





Technical Specifications

Width (base)	550 mm
Depth	880 mm
Height	1770 mm
Safe working load (SWL)	135 kg
Product Code	100060220

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com