


AB2210 Armboard with rotation clamp

- one-handed adjustment armboard with rotation clamp
- versatile adjustment options; height, horizontal and vertical rotation, wide range of angle adjustments
- antimicrobial adjustment handle
- length: 620/ 470 mm
- width: 170 mm
- rod: 150 mm
- Cushion: AB3010 Premium VEF



Technical Specifications

Product code	AB2210 (EU), AB2211 (US), AB2212 (UK)

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com