

Arm and hand table rectangular AB5010

- light and user-friendly arm and hand surgery table with two Easy clamps CL2010
- table top enables fluoroscopy, X-ray imaging area in the middle of table 350 mm
- height adjustable
- rod: 135 mm
- includes mattress





TECHNICAL SPECIFICATIONS

Length	650 mm
Width	400 mm
Product code	AB5010 (EU), AB5011 (US), AB5012 (UK)

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com