

Arm and hand table, curved AB5100

- versatile arm and hand operating table with height-adjustable supporting leg (AB5101) and two Easy clamps CL2010
- provides support for surgeon on each side and allows up to 3 people working around the table
- table top enables fluoroscopy, X-ray imaging area in the middle of table 300 mm
- height adjustable
- rod: 135 mm
- includes mattress





TECHNICAL SPECIFICATIONS

Length	800 mm
Width	600 m
Product code	AB5100 (EU), AB5111 (US), AB5112 (UK)

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com