

## Arm rest 125

Padded arm board 125 includes a fixing clamp and straps 110. Adjustable height, gas spring assisted angle adjustment. Rotates 180 degrees. Suitable for EUR rails.



### Technical Specifications

Width	640 mm
Height	205 mm
Thickness	60 mm
Max. patient weight	325 kg
Safe working load (SWL)	19.5 kg
Trendelenburg	-36 °
Reverse Trendelenburg	+26 °
Weight	6.2 kg

Product Code

110125

## Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)