

## Arm rests 151

Arm rest pair 151 for ophthalmology is designed to support surgeon's arms during operation. Arm board turns 360 degrees. Height adjustment with joint screw on the side of the aluminium clamp. Fits all Merivaara operating tables. Rod diameter 16 mm. Including clamps for fixing. Suitable for EUR, US and UK rails.



### Technical Specifications

|                         |              |
|-------------------------|--------------|
| Width                   | 300 mm       |
| Length                  | 85 mm        |
| Height adjustment       | 120 - 320 mm |
| Thickness               | 40 mm q      |
| Safe working load (SWL) | 10 kg        |
| Weight                  | 2.4 kg       |
| Product Code            | 100000151    |

## Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)