

Armboard long rod AB2130

- long rod model
- one-handed adjustment armboard with Easy clamp
- wide range of height and angle adjustment with ball joint
- antimicrobial adjustment handle
- rod: 280 mm
- includes VEF mattress



TECHNICAL SPECIFICATIONS

Length	620 mm
Width	170 mm
Product code	AB2130 (EU), AB2131 (US), AB2132 (UK)

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com