

Hand fixation device 10393

- Weinberger-type hand fixation device
- For hand and finger traction
- Can be used only for hand operating table 10390
- Translucent and sterilizable.
- Size: 180 mm x 179 mm
- Weight: 0,3 kg



Technical Specifications

| | |
|--------------|------------------|
| | |
| Weight | 0,3 kg |
| Product code | 10393, 100010393 |

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com