

## Hand operating table 10390

Special shaped radiolucent hand operating table with one support leg. Its curved design allows for the surgeon to get closer to the operating area without the need to lean forward. Suitable for EUR, US and UK rails. Includes mattress and clamps. Option to use traction device 10392 and hand fixation device 10393 (Weinberger-type).



### Order no.

110390	Hand operating table; Options: EU, US, UK rails
100 010 392	Traction device
100 010 393	Hand fixation device (Weinberger-type)

### Technical Specifications

Width	850 mm
Height	520 mm
Height adjustment	700 - 1150 mm
Thickness	30 mm
Max. patient weight	200 kg
Safe working load (SWL)	15 kg
Weight	8.8 kg
Product Code	110390, 100010390

## Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)