


Multi-task armboard MV-11802

Arm rest provides arm support for prone, neuro, lateral and chair positions. A wide range of position flexibility, turns 360 degrees and is height adjustable. The premium durasoft pad ensures patient comfort. Universal clamp 10300 is needed for fixing.

Order no.

| | |
|-------------|-----------------------------------|
| 110802 | Multi-task armboard MV-11802 with |
| 100 011 802 | Lateral arm board MV-11802 |
| 100 010 300 | Universal clamp 10300 |



Technical Specifications

| | |
|-------------------------|--------------|
| Width | 510 mm |
| Length | 150 mm |
| Height adjustment | 170 - 480 mm |
| Thickness | 30 mm |
| Max. patient weight | 225 kg |
| Safe working load (SWL) | 16 kg |
| Weight | 2.2 kg |
| Lateral tilt | 45 - 45 ° |

Product Code

110802, 100011802

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com