

## Arm and hand table, curved AB5100

- versatile arm and hand operating table with height-adjustable supporting leg (AB5101) and two Easy clamps CL2010
- provides support for surgeon on each side and allows up to 3 people working around the table
- table top enables fluoroscopy, X-ray imaging area in the middle of table 300 mm
- height adjustable
- rod: 135 mm
- includes mattress





## TECHNICAL SPECIFICATIONS

Length	800 mm
Width	600 m
Product code	AB5100 (EU), AB5111 (US), AB5112 (UK)

### Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)