

## Armboard long rod AB2130

- long rod model
- one-handed adjustment armboard with Easy clamp
- wide range of height and angle adjustment with ball joint
- antimicrobial adjustment handle
- rod: 280 mm
- includes VEF mattress



### TECHNICAL SPECIFICATIONS

|              |                                       |
|--------------|---------------------------------------|
|              |                                       |
| Length       | 620 mm                                |
| Width        | 170 mm                                |
| Product code | AB2130 (EU), AB2131 (US), AB2132 (UK) |

### Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)