

## Easy clamp Ø 18 mm

- light and user-friendly aluminium lever clamp
- antimicrobial tightening levers
- for attaching rods with a diameter of 16-18 mm



### TECHNICAL SPECIFICATIONS

Product code	CL2010 (EU), CL2011 (US), CL2012 (UK)

### Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)