

Easy clamp Ø 20 mm

- light and user-friendly aluminium lever clamp
- antimicrobial tightening levers
- for attaching rods with a diameter of 20 mm



TECHNICAL SPECIFICATIONS

| | |
|--------------|---------------------------------------|
| | |
| Product code | CL2020 (EU), CL2021 (US), CL2022 (UK) |

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com