

## Easy clamp Ø20 mm rotation

- light and user-friendly aluminium rotation lever clamp
- rotation feature allows versatile functionality and adjustability
- antimicrobial tightening levers
- for attaching rods with a diameter of 20 mm



### TECHNICAL SPECIFICATIONS

Product code	CL2120 (EU), CL2121 (US), CL2122 (UK)

### Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)