

## Easy clamp flat bar

- light and user-friendly aluminium lever clamp
- antimicrobial tightening levers
- for attaching flat rods with dimensions of 3-7 mm x 25 mm



## TECHNICAL SPECIFICATIONS

Product code	CL2210 (EU), CL2211 (US), CL2212 (UK)

### Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)