

# DORO Lucent Head Rest System

- Optimized for iCT imaging procedures
- No lateral obstruction due to no overhang at the Base Unit
- Multiple easy ways of mounting the Base Unit
- Long reach into gantry
- Design with a clear focus on usability, cleanability and radiolucency.

## Includes the items below:

110 102 100	Storage Case (Head rest comes in Base Unit
110 103 100	Locking Transitional Member
110 102 600	Transitional Member
100 303 300	Swivel Adaptor Radiolucent
100 303 400	Skull Clamp Radiolucent
110 600 300	Disposable Skull Pins, Adult



## Technical specification

Product code	110102000
--------------	-----------

## Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)

