

ENT-surgery head rest 18143

Gas assisted narrow head rest with contour moulding for comfort, vertical up and down manipulation.

Fits all Merivaara operating tables with "Beach Chair" back rest. For Promerix and Practico the fixing adapter 60146 is needed.

Order no.

110143 Head rest 18143 with adapter 60146 100
018 143 Head rest 18143 100 060 146 Adapter
for Promerix/Practico



Technical Specifications

Width	242 mm
Height	245 mm
Thickness	40 mm

Max. patient weight	325 kg
Safe working load (SWL)	14 kg
Trendelenburg	-45 °
Reverse Trendelenburg	+25 °
Weight	3.4 kg
Product Code	110143, 100018143

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com