

Anaesthesia screen 111

Adjustable height and angle. Rotates 360 degrees. Made of stainless steel. Universal clamp 10300 is needed for fixing. Rod diameter Ø 19 mm

Order no.

110111	Anaesthesia screen kit, incl. anaesthesia screen 111, 2 x extension arms 116 and 10300 clamp (for EUR, US and UK rails)
100 000 111	Anaesthesia screen 111
100 000 116	Extension arm for anaesthesia screen 116



Technical Specifications

Length	590 mm
Width (base)	20 mm
Weight	1 kg
Product Code	110111, 100000111

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com