

## K-Tek 250 stirrup pair 12814

The Kyra 250 is a junior stirrup designed for patients aged 6-12 years old but also effective in smaller/lighter adults weighing up to 250 lbs/113 kg.

- More working space with wider boot rotational range
- Lightweight and Strong with Kyra's Core-Grip™ technology

### Order information

100012814	K-Tek 250 stirrup pair
100012808	K-Tek™ stirrup cart
100012811	K-Tek™ stirrup boot pad kit
100 011 043	Rail clamp, EUR and UK
100 011 043U	Rail clamp, US



## Technical Specifications

Max. patient weight	113 kg
Weight	4.6 kg
Movement angle	+ -25
Lithotomy Range	+85 - -55 °
Abduction range	+25 - -9 °
Product Code	100012814

### Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)