

K-Tek 350 stirrup pair 12806

The K-Tek Comfort Stirrups for patient positioning for a wide range of gynecologic, urologic and laparoscopic procedures. The Stirrups increase surgeons' operating room performance by permitting simple adjustment of lithotomy and abduction intra-operatively. The K-tek Comfort 350 stirrup is designed for patients with a BMI of less than 29 with a large flexible lateral fin and can support patients up to 159 kg weight capacity.

- Designed to relieve pressure with soft padding that encapsulates the foot.
- K-TEK™ silicone closure straps making cleaning easy.
- Flexible finned boots accommodates various leg sizes.
- Superior access with lithotomy range of motion from -55° to $+85^{\circ}$ aids in relieving clinician shoulder/arm strain.



Order information

100012806	K-Tek Comfort 350 stirrups
100012808	K-Tek™ stirrup cart
100012811	K-Tek™ Comfort™ 350 stirrup additional boot pad kit (one piece)
100 011 043	Rail clamp, EUR and UK
100 011 043U	Rail clamp, US



Technical Specifications

Max. patient weight	159 kg
Lithotomy Range	+85 - -55 °
Abduction range	+35 - -9 °
Product Code	100012806

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com