

# Kyra Blue 400 Stirrups with Secure-Lok

The K-tek 400 stirrup is designed for bariatric patients.

- Non-finned boot
- Can be used with robotic surgery
- Easy to clean durable silicone straps
- Includes more padding to reducing pressure
- Increased range of motion



## Order information

100 012 815	K-Tek 400 stirrup pair
100012808	K-Tek™ stirrup cart
100012811	K-Tek™ stirrup boot pad kit
100 011 043	Rail clamp, EUR and UK
100 011 043U	Rail clamp, US

## Technical Specifications

Max. patient weight	181 kg
Product code	100012815

## Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)