

## Leg rest, pair 20118

Leg rests 20118 to be used with Lithotomy or as an over arm rest when patient is in lateral position. Goepel type leg rests are padded and with velcro straps. Wide trajectory thanks to its ball joint. Stainless steel. Including universal clamps 10300 (2 pcs).

### Order codes:

110118	Leg rests with straps 20118, pair, incl. 2 x 10300 clamps for EUR, US and UK rails
100 020 118	Leg rest 20118 (pair)
100 010 300	Universal clamp (1 pc)



### Technical Specifications

Length	320 mm
Height	650 mm
Max. patient weight	225 kg
Safe working load (SWL)	29 kg
Weight	5,1 kg
Product Code	110118, 100020118

## Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)