

Femoral counter traction post 519512

Femoral counter traction post for lateral position is to be used with orthopedic extension device Flexion. Femoral counter traction post is meant for femur treatment, medullary nailing in lateral position. The post is equipped with a height adjustment and removable curved cross bar with conical pads.





Technical Specifications

| | |
|-----------------------------|-----------|
| | |
| Безопасная рабочая нагрузка | 135 kg |
| код продукта | 100519512 |

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com