


## Knee support set

Knee support with carbon fibre bar and mattress is meant for knee surgery. It is ideal to be used with orthopedic extension devices for tibia nailing. Radial universal clamp 100010300 for fixing included.

- Diameter Ø 18 mm
- Size 290 mm, Ø 120 mm
- Bar length 581 mm
- Mattress thickness 35 mm
- Order No. A43134200



### Technical Specifications

Thickness	35 mm
Product Code	A43134200

## **Contact us**

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)