

## Short support arm 19313

Short support arm 19313 for supporting healthy leg with stirrup or leg support, e.g. 20118. Clamp 19135 to be ordered separately.



### Technical Specifications

Length	50 mm
Height	300 mm
Max. patient weight	225 kg
Safe working load (SWL)	58 kg
Weight	1.5 kg
Product Code	100019313

## Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)