

# Trenguard Hybrid

## For Trendelenburg patient restraint

TrenGuard is an advanced patient positioning solutions for gynecological and urological procedures and robotic surgery.

### Key benefits

- No sliding – patented “Speed Bump” bolster keeps patients of all sizes from sliding
- No shoulder braces
- Soft lateral stabilizing pillows control body mass shift
- No chest straps
- Promotes normal ventilation
- Efficient and compatible
- Accommodates under-patient devices
- Simplifies patient repositioning
- Less waste
- Re-usable Lateral stabilizing pillows mean less medical waste
- Lower cost per procedure
- System easily accommodates patients up to 250 kg



**450 Hybrid starter package Standard**

Including:

- Patient Support Frame
- Clamps for EUR rails
- Single-use soft head stabilizing pillow and Speed Bump bolster
- Lateral Stabilizing Pillows
- Product code A43824200

## Reusable package

Including:

- One reusable soft lateral stabilizing pillow control body mass shift during transition to Trendelenburg
- Order No. 55251

## Reusable package for bariatric patients

Including:

- One reusable soft lateral stabilizing pillow to control body mass shift during transition to Trendelenburg
- Order No. 56251

## Single-use wedge package

For bariatric patients and those without discernible cervical notch.

Including:

- 12 single-use integrated wedge and head stabilizing pillows
- Order No. 56401

## Single-use package

Including:

- 12 single-use Speed Bump bolsters which prevent patient sliding
- 12 single-use soft head stabilizing pillows which control side-to-side rocking of the head
- Order No. 55201

## Single-use package for bariatric patients

Including:

- 12 single-use Speed Bump bolsters which prevent patient sliding
- 12 single-use soft head stabilizing pillows which control side-to-side rocking of the head
- Order No. 56201

## Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)