

## Knee support SS2040

- support roll for knee and leg surgery treatment
- suitable also with orthopaedic traction device
- height adjustment range 600 mm
- pad diameter 100 mm, length 380 mm
- Recommended clamp: CL2110



### TECHNICAL SPECIFICATIONS

Product code	SS2040

### Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)