

## Multi lateral support, wide SS1050

- multi-adjustable lateral support for optimum positioning
- padding (50 mm) designed to prevent pressure and nerve injuries
- antimicrobial tightening levers
- height adjustment range 200 mm
- pad size 200x200 mm
- Recommended clamps: CL2010/ CL2110



### TECHNICAL SPECIFICATIONS

Product code	SS1050

### Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)