

Floating Finn Stirrups

The new Merivaara Floating Finn lithotomy stirrups are designed for effortless and optimal lithotomy positioning. They are available in two sizes: Floating Finn for adults with a safe working load (SWL) of 180 kg and Floating Finn Plus for bariatric patients with a SWL of 250 kg.

- Non-Velcro silicone straps support improved infection control.
- Flexible boot accommodates larger patient calves.
- Pressure-relieving padding reduces the risk of pressure injuries and compartment syndrome.
- Clearly visible scale simplifies symmetrical alignment.



Order information

100012850	Merivaara Floating Finn Stirrups, including Boot Pads - Silicon Straps
100012851	Merivaara Floating Finn PLUS Stirrups, including Boot Pads - Silicon Straps
100012854	Merivaara Floating Finn Cart
100012853	Merivaara Floating Finn Stirrup Replacement Pads
100012852	Merivaara Floating Finn Plus Stirrup Replacement Pads



Technical Specifications

	Floating Finn Stirrups	Floating Finn PLUS
Safe working loag (SWL)	180 kg	250 kg
Product Code	100012850	100012851

Связаться с нами

Merivaara Corp.

Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com